


Note: Part numbers and standards' values may change. Consult LECO for the latest information.



**GENERAL COMBUSTION SUPPLIES**

621-331	Ceramic Crucible—Large (5 pkg)
529-048	Ceramic Crucible Cover—Large (5 pkg)
502-628	Crucible Foil Liner (1200/pkg)
761-929	Crucible Tongs
778-898	Crucible Cover Forceps
614-442-110	Crucible Turntable—Large (Kit)
776-978	Glass Scoop
777-046	Heat Resistant Glove

**OPTIONS**

621-010-140	Balance Calibration Kit (100 g) Includes: 100 g weight and 4-inch forceps
619-592-407	Air Regulator
619-592-406	Inert Gas Regulator
619-592-403	Oxygen Regulator
621-225G	External Blower Kit (115V)
621-361	~5 g Quartz Test Rod

**REFERENCE MATERIALS****COAL PROX-PLUS (50 grams)**

Part No.	Desc.	% Volatile	% Ash
502-680	Low	~17.0	~11.0
502-681	Medium	~32.0	~9.0
502-682	High	~41.0	~6.0

**COKE PROX-PLUS (50 grams)**

Part No.	Desc.	% Volatile	% Ash
502-684	Pet Coke Prox Plus	~11.7	~0.4
502-683	Met Coke Prox Plus	~1.0	~10.0

**CALCIUM OXALATE (50 grams)**

Part No.	% Weight Loss		
	@ 200°C	@ 400°C	@ 850°C
502-091	12	19	30

*NOTE: Part numbers and standards' values may change.  
Consult LECO for latest information.*